

Table Top Autoclave

Model : **ESL-25D , ESL-50D , ESL80D**



Features

- Full automation
- LED Display
- Safety device
- Self-test function
- Exhaustive design's gasket



MODEL		ESL-25D	ESL-50D	ESL-80D
Dimension (W x D x H) (mm)	(In)	535 x 620 x 430	620 x 710 x 520	620 x 1000 x 1290
	(Out)	Ø260 x 460mm	Ø350 x 510mm	Ø1350 x 800mm
Capacity		25 Liter	50 Liter	80 Liter
Heater		2KW	2KW	2.5KW
Weight		48kg	65kg	100kg
Tray size		190 x 400 x 25mm	240 x 480 x 30mm	240 x 770 x 30mm
Tray material		Aluminum		
Reservoir capacity		4.5 Liter	4.5 Liter	8 Liter
Material	(In)	Stainless steel		
	(Out)	Steel plate with powder coating		
Safety device		Over temp. protector, Safety Valve		
Temp.range		121°C		
Pressure range		0-2.1kgf/cm ²		
Sterilizer Time		60Min(Max)		
Control system		8 Bit Microprocessor		
Power		220V, 50/60Hz		

