

Shaking Incubator

Model **ESL-8480SF**
ESL-8480S

Labware



Main Feature

- Noiseless & Brushless DC Motor
- Durable & Beltless Drum type
- Stable operation in low speed of 10RPM
- Digital PID Controller
- Reliable safety device
- Simple-to-use
- Compatible for various capacity of flask & tube
- Light Bank ("-L" Model)

Shaking Incubator

Model	ESL-8480SF	ESL-8480S
Temp. Range	15°C to 60°C	Ambient +5°C to 60°C
Weight	81Kg	66Kg
Power	220V, 4.5A , 60Hz	220V, 1.8A, 60Hz
Temp. Controller	Digital P.I.D. controller	
Speed Range	10 - 350 RPM	
Display	LED Display	
Drive system	Beltless direct drive	
Motor	Brushless D.C.Motor K-type	
Circulation Fan	3.5" convection fan x 2EA	
Dimension (IN)	510 x 530 x 400(H)mm	
(OUT)	515 x 730 x 720(H)mm	
Lid	Acrylic Blue	
Safety Device	Over temp.Protector, Door detector	
Standard plate	250ml x 25EA	

Light Bank Model

Model	ESL-8480SFL
Light Bank	40W CL x 1 Unit
Light Control	24Hr. Photoperiod program

Plate of Flask

100ml Flask	36 EA
250ml Flask	25 EA
500ml Flask	16 EA
1000ml Flask	9 EA
2000ml Flask	4 EA
Spring Rack	450 x 450mm

